

FLOSSERPIK POWER SPA®



FLOSSERPIK® POWER SPA IS A EFFECTIVE WATER FLOSSER THAT CLEANS THOSE 'HARD TO REACH AREAS'.

Unique Flosserpiik® Power Spa - The ultimate tool for superior oral health.

Effective * Waterproof * Convenient * Compact * Cordless * Rechargeable * Portable

Versatile

Flosserpiik® Power Spa combines a selection of specialist attachments with various cleaning modes and pressure settings.

Portable

Flosserpiik® Power Spa has a slim and lightweight design. It has a long-lasting battery life, so only needs charging once every 1-2 weeks.

Pulsating

Flosserpiik® Power Spa is equipped with a pulsating engine, which not only reduces water consumption, but also provides a more thorough clean.

Waterproof

Flosserpiik® Power Spa is fully waterproof for up to 1m. Ideal for use in the shower!

Easy to clean

The Flosserpiik® Power Spa tank opens from the bottom, allowing easy cleaning.

Mouthwash compatible

Compatible to use with **Perio Plus** Mouthwash (excluding mouthwash containing iodine, hydrogen peroxide or tea tree oil).



- Water flossers remove 99% of plaque.
- Cleans Periodontal pockets of 6mm effectively.
- Proven to significantly improve gingival health.
- Easier, quicker, and more effective than flossing.
- Water jets up to 9m.
- The water tank can be fully opened and cleaned thoroughly.
- Waterproof upto 1m. Ideal for in-shower use.
- Compatible with **Perio Plus** mouthwash to boost freshness and/or deliver oral medication.

PERIOPLUS+



4 EXPERT ATTACHMENTS

Jet tip

Pointed attachment that emits a powerful jet of water. Ideal for interdental and gingival cleaning, as well as water flossing around implants and bridges.

Tongue scraper


Flat, wide head for water-powered tongue cleaning, boosting bacterial balance.

Subgingival tip

Tiny-tipped precision attachment specifically for cleaning below the gum line and deep into periodontal pockets. Use only on Sensitive mode, and as advised by your dentist.

Orthodontic tip

Water-assisted toothbrushing with signature UltraSoft Konex HD® bristles. This dual-cleaning tip flosses while you brush. Perfect for cleaning around brackets and wires of braces.



3 EASY CLEANING MODES

Active mode

- For deep cleaning around teeth and gums
- Ideal for effective teeth brushing and cleaning braces (with the orthodontic tip) or cleaning interdentially and around appliances with the jet tip
- Can be used with any of the attachments
- 2 minutes of continuous pressurised water spraying

Sensitive mode

- Especially suited to clean tender areas or for cleaning periodontal pockets with the subgingival tip.
- 2 minutes of continuous low pressure water spraying

Pulse mode

- For powerful cleaning between teeth
- Ideal for deep cleaning with the jet tip or orthodontic tip
- 3 sprays per second for 2 minutes

